

Test Report

No. BOYNPJIR828897L1a

Sample Description

Medical Mask (Non-Sterile)

Applicant

Shandong Zhushi Pharmaceutical Group Co.,
Ltd..


Pony Testing International Group
www.ponytest.com



Test Results

No. BOYNPJIR828897L1a

Page 1 of 2

Sample Description	Medical Mask (Non-Sterile)	Sample Specification	---
Applicant	Shandong Zhushi Pharmaceutical Group Co., Ltd.	Trade Mark	---
Received Date	2020-03-27	Manufacturing Date or Lot No.	2020.03.02
Test Date	2020-03-27~2020-04-09	Sample Grade	---
Sample Quantity	50pcs	Test Type	Commissioning Test
Sample Status	Flat Face Mask	Test Environment	To meet the requirements
Test Items	See the next page		
Reference Methods	See the next page		
Main Instruments	Electric Heating Constant Temperature Incubator、 Bacterial filtration efficiency detector etc.		
Note	Limit on: EN 14683: 2019 Sample sources: Self produced This test report is to replace the test report No. BOYNPJIR828897L1 The No. BOYNPJIR828897L1 test report is invalid and of no legal effect. All related information should be referred to the new test report. 2020-04-14		
	Edited by	Bianroujia	
	Checked by	swjide Shero	
	Approved by	Dai Qiong	
	Issued Date	2020-04-14	

Test Results

No. BOYNPJIR828897L1a

Page 2 of 2

Sample Description and Number	Test Items	Limit			Test Results		Reference Methods
		Type I	Type II	Type II R			
R828897L1 Medical Mask (Non-Sterile)	Bacterial filtration efficiency(BFE), %	≥ 95	≥ 98	≥ 98	No.1	99.1	EN 14683: 2019
					No.2	99.2	
					No.3	99.0	
					No.4	99.3	
					No.5	99.2	
	Differential pressure, Pa/cm ²	<40	<40	<60	55.5		EN 14683: 2019
	Splash resistance pressure, kPa	/	/	≥ 16.0	> 16.0		EN 14683: 2019
	Microbial cleanliness, cfu/g	≤ 30	≤ 30	≤ 30	<1		EN 14683: 2019

Sample Number and Photo:



Pony authenticate the photo on original report only

—End of Report—

